

Technical Data Sheet

Cryogenic tubes pre-filled with Cytodelics Stabiliser



1. Shipping

Way of transportation	Minimal conditions to meet	Optimal conditions to meet
Road	<ul style="list-style-type: none">• Shipping on icepacks• Packing in plastic bags• Temperature under 25°C	<ul style="list-style-type: none">• Shipping on icepacks• Packing in plastic bags• Temperature under 15°C
Air	<ul style="list-style-type: none">• Shipping on icepacks• Packing in plastic bags• Upright position• Temperature under 25°C	<ul style="list-style-type: none">• Shipping on icepacks• Packing in plastic bags• Upright position• Temperature under 15°C

2. Storing

Maximal stability of the "Cytodelics **Stabiliser**" pre-filled in a cryogenic tube can be reached by refrigeration in the range **from +5°C to +8°C**.

Cryogenic tubes pre-filled with the Cytodelics Stabiliser can be stored at room temperature (20-25°C) for up to 12 weeks.



3. Dimensions

Tube volume	Total tube		Tube without the cap	Space for sticker		Space for writing		Cap height in thread	Rim of the cap to bottom
	Diameter [mm]	height [mm]		height [mm]	height [mm]	Height [mm]	Width [mm]		
1.2 ml	12.5	40.5	34.5	19	17	19.5	13.5	6	27
2.0 ml	12.5	48	41.5	27	21	19.5	13.5	6	34.5
3.0 ml	12.5	69.5	63.5	49	36	19.5	13.5	6	56
4.0 ml	12.5	75	69	54	40	19.5	13.5	6	61.5
5.0 ml	12.5	89.5	83.5	69	55	19.5	13.5	6	76

